# IMPORTANT SAFETY INSTRUCTIONS

## DO NOT DESTROY

**DANGER – To reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions:**

1. Read all instructions.
2. **ALWAYS UNPLUG THE UNIT FROM THE ELECTRICAL OUTLET IMMEDIATELY AFTER USING.**
3. Do not use while sleeping. Avoid falling asleep while using this product.
4. Do not use on an infant or an animal.
5. This pad is not to be used on or by any individual with certain disabilities, a sleeping or unconscious person, or a person with poor blood circulation unless carefully attended.
6. Do not use on areas of insensitive skin.
7. Burns can occur regardless of control setting; check skin under pad frequently.
8. Never use pad without removable cover in place.
9. Do not use where oxygen is being administered.
10. Do not use pins or other metallic means to fasten this pad in place.
11. Do not sit on or crush pad – avoid making sharp folds in pad.
12. Never pull this pad by the line cord.
13. Do not use the cord as a handle.
14. Carefully examine inner cover and line cord before each use. Discard the pad if inner covering or line cord shows any sign of deterioration.
15. Do not tamper with this pad in any way. There are no user-serviceable parts. If, for any reason, this pad does not function satisfactorily, return it to authorized service center at address given.

## HEATING PAD RESTRICTIONS

For additional warnings, please see Important Safety Instructions as well as Cautions in this booklet.

Not for use by children

Do not sleep while using heating pad

Do not use heating pad bunched up

Do not use under covers

Do not sit on heating pad

Do not sharply fold or crease heating pad

16. Do not use this pad with liniment, salve or ointment preparations that contain heat-producing ingredients, as skin burns could result. Use wherever hot applications are desirable for personal comfort and wherever recommended by your physician or dentist for the relief of pain.
17. Loop cord loosely when storing; tight wrapping may damage cord and internal parts.
18. Keep the pad away from heated surfaces.
19. Place pad on top and not under the part of the body needing heat.
20. Do not put this pad in the microwave.
21. Do not leave pad unattended, especially if children are present.
22. Do not use pad underneath bedcovers and other types of blankets. Do not place covers, such as sheets and blankets, over pad. Only the cloth cover provided with the pad should be used.
23. Do not use pad on sensitive skin or when taking medication that can cause sensitive skin.
24. Do not exceed 30 minutes of heat treatment or 10 minutes of heat and massage treatment at one time.
25. Do not place the unit into a washing machine. Spot clean only, using a soft, damp cloth. Do not use any household cleaners.
26. When not in use, unplug the heating pad and store in a dry location safely out of reach of children.

## FCC MARKINGS

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.

• Consult the dealer or an experienced radio/TV technician for help.  
**Caution:** Any changes or modifications to this unit not expressly approved by Conair LLC, the manufacturer of this product and the party responsible for compliance, could void the user's authority to operate the equipment

## SUPPLIER'S DECLARATION OF CONFORMITY

**Brand:** CONAIR  
**Model No.:** HP08FX  
**Description:** MASSAGING HEATING PAD  
**Responsible Party:** Conair LLC,  
 1 Cummings Point Rd.  
 Stamford, CT 06902  
 (203) 351-9000

**Standards:** FCC Part 15 Subpart B, ANSI C63.4-2014  
 This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesirable operation.

**DATE OF ISSUE:** FEBRUARY 22, 2021

# SAVE THESE INSTRUCTIONS

## Cautions:

- **This product does not have an automatic-off for the heat function. This product has a 20 minutes automatic-off feature for massage.**
- This heating pad should not be used in the following cases unless directed by a physician: if you have diabetes, tuberculosis, benign or malignant tumors, phlebitis or thrombosis, hemorrhages, inflammations of the skin, undetermined calf pain or joint pain, or poor circulation. Consult a physician before using if you have any kind of fungal infection, swelling, fracture or persistent pain.
- **Consult your physician if you have any questions regarding therapeutic use of this heating pad.**
- Do not use on any unexplained pain or swollen muscles or following a serious muscular injury before consulting your physician.
- **Keep away from children.**
- Discontinue use of product if pain increases or swelling occurs after application.
- Do not use while bathing or showering. Use only in a dry environment.
- Consult your doctor before using in cases of illness or medical condition, including but not limited to pacemaker use, pregnancy or cancer.
- If you experience pain in a muscle or joint for a prolonged period of time, discontinue use and consult your doctor. Persistent pain could be a symptom of a more serious condition.
- Never use any heating products on open wounds, discolored areas, or any area of the body that is swollen, burned, bruised or inflamed or where skin eruptions, ulcerated sores or varicose veins are present.

- This product should never be used by any individual suffering from any ailments that would limit or altogether inhibit the user's capacity to feel or have any sensation in the body part that this product is being used for.
- Use of this product should be pleasant and comfortable. Should pain or discomfort result, discontinue use and consult your physician.
- Ensure that your hands are dry when operating the controls or removing the plug.
- Do not fall asleep while using this product.
- For indoor use only.
- **DO NOT EXCEED 30 MINUTES OF HEAT TREATMENT OR 10 MINUTES OF HEAT AND MASSAGE TREATMENT AT ONE TIME.**

## OPERATING INSTRUCTIONS

To operate the ConairComfort™ Massaging Heating Pad, first plug the power jack into the hand control then plug the cord into a standard 120V electrical outlet.

1. Press the Heat button (bottom button) to turn the heat on or off.
2. Press the Massage button (top button) repeatedly to adjust the level of massage to High or Low, and to turn the massage feature off.

The following light indicators show which function(s) is/are being used:

<b>Red indicator</b>	<b>Heat only</b>
<b>Red + Green indicator</b>	<b>Heat + Low Massage</b>
<b>Red + Orange indicator</b>	<b>Heat + High Massage</b>
<b>Green only</b>	<b>Low Massage Only</b>
<b>Orange only</b>	<b>High Massage Only</b>

For maximum benefit, use should not exceed 30 minutes of heat, or 10 minutes of heat and massage at one time. The Massage function will automatically shut off after 20 minutes of use.

For optimal use, place the side of the heating pad with the "Heated Side" label closest to your body. If the heating pad becomes too hot, flip the pad over and use the opposite side

## LIMITED ONE-YEAR WARRANTY (U.S AND CANADA ONLY)

Conair will repair or replace (at our option) your unit free of charge for 12 months from the date of purchase if the appliance is defective in workmanship or materials.

To obtain service under this warranty, return the defective product to the service center listed, together with your purchase receipt. California residents need only provide proof of purchase and should call 1-800-3-CONAIR for shipping instructions. In the absence of a purchase receipt, the warranty period shall be 12 months from the date of manufacture.

**ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 12-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY.** Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

**IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER.** Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you.

This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

Please register this product at [www.conair.com/registration](http://www.conair.com/registration)

**Service Center**  
 Conair LLC  
 Service Department  
 7475 North Glen Harbor Blvd.  
 Glendale, AZ 85307

stock.image

Visit us on the web at [www.conair.com](http://www.conair.com)

©2021 Conair LLC  
 Glendale, AZ 85307

HP08FX  
 21PA079517

# CONAIRcomfort™ massaging heating pad

## instruction booklet

for Model HP08FX

For your safety and continued enjoyment of this product, always read the instruction book carefully before using.

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